


# OM-Mount Photomicro Adapter H + PM-D35-S



## OM-Mount Photomicro Adapter H + PM-D35-S

- This adapter is mainly used to mount OM System SLR cameras for photography in combination with Olympus special photomicrographic system.
- When mounting the OM camera body on the Photomicrographic System PM-10, automatic or manual version, or the Photomacrographic System PMT-35A, the Adapter H is needed together with the Camera Adapter PM-D35-S.
- When photographing with the PM-10 Series Photomicrographic System the shutter of the OM camera body should be kept open at the "B" Bulb setting and the vibration-free shutter built into the PM-10 is used to make the exposure. This will assure sharp pictures free of camera shake even at high magnifications of 40x or more.



OM-Mount Photomicro Adapter H Stock No. 105-037  
 35mm SLR Camera Adapter PM-D35S Stock No. 105-037

**OLYMPUS**

# OM SYSTEM OM-Mount Photomicro Adapter H + PM-D35-S

## Meter Coupling

Sets the exposure meter built into the camera body to the exposure value necessary for photomicrography.

## Camera Mounting Ring


Bayonet mount connection to OM body as with ordinary lenses.

## Adapter Release Button

Press the button and rotate the Adapter H 70° counterclockwise to remove. Do not press button when mounting.

## Clamping Screw


After inserting the Adapter H, tighten this screw to secure.



## 35mm SLR Camera Adapter PM-D35-S

Connect to the Photomicrographic System PM-10 Series Exposure Body or Macrophotographic System PMT-35A.

## Typical use of OM-Mount Photomicro Adapter H and PM-D35-S



## Notes

- With the BH-2 Series, precise focusing can be done through the binocular tube, so it is not necessary to use the Photomicrographic Equipment Focusing Telescope or the SLR viewfinder.
- With microscopes other than the BH-2 Series (the CH, CK2, various stereoscopic microscopes, etc.), focusing should be done with the Focusing Telescope of the Photomicrographic Equipment.
- When photographing with the PMT-35A, focusing is accomplished through the SLR camera viewfinder. The Focusing Screen 1-12 (P. C-001), which is specially designed for photomicrography, should be used. For sharper focusing and a more comfortable viewing position, use of the Varimagni Finder (P. C-003) together with the Eyecoupler (P. C-002) is recommended.